

Scent trails: tracking HS production with the MicroMap

Author: Cyrille C. Thinnés, University of Galway

This case study presents a worked example of MicroMap applications. It introduces key concepts in microbiome metabolism and constraint-based modelling, serving as both an educational resource and an accessible entry point to the COBRA Toolbox [1]. You can find an [accompanying video walkthrough on YouTube](#). We recommend first completing the '*MicroMap exploration with CellDesigner*' tutorial.

Introduction

The most abundant gases contributing to flatulence are odourless, such as hydrogen (H₂), carbon dioxide (CO₂), and methane (CH₄) [2]. The odoriferous onus pertains to trace gases (<1% of total volume), notably sulphur-containing compounds, of which hydrogen sulphide (HS) has been identified to be the most putrid perpetrator [3]. HS is produced from microbiome metabolism and is linked to health and disease [4]. We shall use the MicroMap to gain mechanistic insights into the metabolism of this vicious volatile.

MicroMap inspection with CellDesigner alone

First, we inspect the MicroMap in CellDesigner to gain preliminary mechanistic insights into microbial HS metabolism. Open the MicroMap and search for HS using Ctrl+F (Find Species). Enter the VMH abbreviation h2s, tick name (not ID), and choose 'begin with' to capture all compartments. Press 'Show All'. CellDesigner will highlight all h2s occurrences with pink boxes. When zooming out, you will notice clusters of HS-associated reactions. To aid inspection, recolour these metabolites: with all selected, click 'Change Color & Shape', choose a bright red, enlarge the size (e.g. 200 × 75), and apply to all occurrences. All HS species are now easily spotted across the map.

Zooming out reveals that HS metabolites cluster in the top right quadrant. Since the MicroMap is organised by subsystem, we can now examine which reactions and pathways they belong to. Zooming into the 'Sulfur Metabolism' region shows two distinct HS nodes. Although some reactions are closely positioned and partially overlapping, they can be easily identified by selecting the reaction arrow and reading its name in the 'Notes' window. Reaction appearance can be edited via the 'Change Color & Shape' button (e.g. line width 5, red), and multiple reactions can be highlighted and edited simultaneously by holding 'Shift' while selecting them. We may also highlight reaction cascades by applying colour codes, e.g., by marking reactions directly connected to HS in red and those further up/downstream in orange (**Figure 1**).

We see two distinct sections (**Figure 1**), linked by a common node of sulphite (SO₃²⁻):

- Organic (left): HS is involved in the interconversion of carbon-containing compounds.
- Inorganic (right): HS is involved in the interconversion of sulphur-containing ions.




Figure 1: HS metabolism inspection using the MicroMap. The HS molecules were enlarged and highlighted in bright red to easily identify their location on the map. Reactions directly involving HS were then highlighted in red, with further up- and down-stream reactions in orange, to highlight the wider HS-involving network. This approach reveals an ‘organic’ reaction group primarily involving carbon-containing reaction partners, and an ‘inorganic’ reaction group primarily involving inorganic sulphur compounds.

Proceeding by inspecting the methionine, cysteine, and serine metabolism subsystems reveals recurring amino acid transformations involving HS, including between cysteine and pyruvate, cysteine and serine, homocysteine and cysteine, and acetylserine and cysteine.

In summary, inspection of the MicroMap alone reveals that HS is concentrated in the sulphur, methionine/ cysteine, and serine metabolism subsystems. Multiple direct and indirect routes exist for HS formation, including amino acid transformations and sulphur anion conversions. Across these, the repeated $HS \leftrightarrow SO_2$ interconversion emerges as a central organising feature, linking organic and inorganic sulphur metabolism.

A quick way to assess whether a microbe partakes in HS metabolism is to load its reconstruction visualisation and highlight the HS molecules, as outlined above. If any highlighted HS metabolite is connected to a reconstruction reaction, the microbe can perform that reaction. For example, when comparing the pan-reconstruction visualisations of *Desulfovibrio desulfuricans* and *Lactobacillus pasteurii* ([downloaded from the Harvard Dataverse](#)), the former shows extensive HS metabolism, whereas the latter shows none (**Figure 2**).




Figure 2: Comparison of HS metabolism in different reconstructions. Reactions involving HS can be identified by using the genome-scale metabolic reconstruction visualisation for a given microbe and highlighting all the contained HS metabolites. If HS is connected to an already coloured-in reaction, it will be a reaction partner.

MicroMap integration with the COBRA Toolbox

First steps

The next steps require a current version of the COBRA Toolbox.¹ For installation instructions, please [follow the installation guide](#). Additional information and user documentation are available on the [COBRA Toolbox website](#).

Download the MicroMap CellDesigner .xml file [from the Harvard Dataverse](#) and place it in your chosen MATLAB working directory. If required, initialise the COBRA Toolbox and set the solver.

```
%First Steps
initCobraToolbox(false); %no need to update if you have the latest version
changeCobraSolver('gurobi'); %set your installed solver
```

Parsing the MicroMap into the COBRA Toolbox

CellDesigner .xml files can be parsed in the COBRA Toolbox with the *transformXML2Map* function. This splits the file into two components: the xml part (not relevant here) and the map part, which will be used in subsequent steps.

```
[xmlMicroMap, mapMicroMap] = transformXML2Map('MicroMap.xml');
```

The two parts can later be recombined into a CellDesigner .xml file using *transformMap2XML*:

```
transformMap2XML(xmlMicroMap, mapMicroMap, 'CopyMicroMap.xml');
```

Modifying the MicroMap with the COBRA Toolbox

The MicroMap uses a coherent colour scheme to help recognise related metabolites. For some applications, however, it may be preferable to maximise contrast by simplifying the appearance. The *unifyMetabolicMapCD* function can be used to set all metabolite colours to white and all reaction colours to grey with a line width of 1 (**Figure 3**).

```
%Decolorise the map  
mapMicroMapUnified = unifyMetabolicMapCD(mapMicroMap);  
transformMap2XML(xmlMicroMap, mapMicroMapUnified, 'UnifiedMicroMap.xml');
```




Figure 3: Decolorising the map. The COBRA Toolbox contains the *unifyMetabolicMapCD* function, which enables us to create a blank canvas from an already coloured map. All metabolite nodes will be set to white, and all reaction lines will be set to grey with line width 1.

The metabolic cartography suite includes several functions to adjust metabolite and reaction appearance. For our use case, we want to colour reactions involving HS, which can be done with *modifyReactionsMetabolites*. While the VMH identifier for HS (h2s[c]) is known, the associated reactions must first be identified. The Metabolic Cartography Tutorial Part 2 in the COBRA Toolbox demonstrates how to extract such a list from a genome-scale reconstruction using *findRxnsFromMets*. Here, however, we are working with large resources spanning many reconstructions. One option is to build a pan-reconstruction from all microbiome reconstructions;

alternatively, the full set of VMH reactions can be queried by loading the database with *loadVMHDatabase*. With this in place, we can generate a list of reaction IDs that include h2s[c] as a reaction partner.

```
%Colouring metabolites and reactions.
%Identify all rxns containing h2s as a reaction partner from the VMH

%Load all reactions from the VMH Database
allVMH = loadVMHDatabase;
reactions = allVMH.reactions;

% Locate columns by header names
headers = reactions(1,:);
colFormula = find(strcmp(headers,'Formula'));
colRxnID = find(strcmp(headers,'ReactionID'));

% Data rows
data = reactions(2:end,:);

% Find reaction formulas containing h2s[c] (case-insensitive)
formulas = string(data(:, colFormula));
mask = contains(lower(formulas), "h2s[c]");

% Collect matching reaction IDs
H2Srxns = unique(string(data(mask, colRxnID)));
```

With the necessary inputs collected, we can now use *modifyReactionsMetabolites*. To define specific colours, we shall [employ html colour codes](#). For using, e.g., a red-orange colour with line width of 10, we enter:

```
% Color in H2S reactions
mapMicroMapH2SRxns = modifyReactionsMetabolites(mapMicroMapUnified, ...
    H2Srxns , 'h2s[c]', 'ORANGERED', 10);
transformMap2XML(xmlMicroMap, mapMicroMapH2SRxns, 'H2SRxnsMicroMap.xml');
```

The resulting map highlights all reactions involving HS in microbiome metabolism, mirroring the manually generated CellDesigner map described above (**Figure 2**). Accordingly, the biochemical interpretation remains the same.

Comparing the reaction content in different reconstructions

In the CellDesigner section above, we inspected pre-generated reconstruction visualisations by manually colouring HS metabolites to check their adjacency to reactions of interest (**Figure 2**). With the COBRA Toolbox, we can perform the same analysis from scratch, starting with the reconstructions and the MicroMap and using *visualizeReconstructionsOnMap*. Download the AGORA2 pan-species reconstructions (.mat format) for *Desulfovibrio desulfuricans* and *Lactobacillus pasteurii* [from the Harvard Dataverse](#) (you can search with wildcards, e.g. *desulfuricans). Place the reconstructions in the same directory (not necessarily the working directory). For visualisation, run:

```
% Visualising reconstructions
```

```
visualizeReconstructionsOnMap('MicroMap.xml', pwd);
```

(assuming both the MicroMap and reconstructions are in your MATLAB working directory – you can also specify exact paths if needed.) Running the code will generate visualisations for each reconstruction: `MicroMap_panDesulfovibrio_desulfuricans.xml` and `MicroMap_panLactobacillus_pasteurii.xml`. Metabolite colours can be modified with `changeMetColor`, and node sizes with `changeNodesArea`. Remember to parse each map into MATLAB before making changes (see above). Thus, we can run:

```
%Parse both maps into Matlab
[xmlDDesu, mapDDesu] =
transformXML2Map('MicroMap_panDesulfovibrio_desulfuricans.xml');
[xmlLPast, mapLPast] =
transformXML2Map('MicroMap_panLactobacillus_pasteurii.xml');

%Change the colour for the h2s metabolites
mapDDesuH2S = changeMetColor(mapDDesu, {'h2s[c]'}, 'LIME');
mapLPastH2S = changeMetColor(mapLPast, {'h2s[c]'}, 'LIME');

%Change the node size for the h2s metabolites
mapDDesuH2S = changeNodesArea(mapDDesuH2S, {'h2s[c]'}, 75, 150);
mapLPastH2S = changeNodesArea(mapLPastH2S, {'h2s[c]'}, 75, 150);

%Reconstitute into CellDesigner maps
transformMap2XML(xmlDDesu, mapDDesuH2S,
'MicroMap_H2SpanDesulfovibrio_desulfuricans.xml');
transformMap2XML(xmlLPast, mapLPastH2S,
'MicroMap_H2SpanLactobacillus_pasteurii.xml');
```

With this approach, it is straightforward to see which reconstructions contain reactions involving HS by checking whether a highlighted HS node is linked to a highlighted reconstruction reaction in CellDesigner. In this example, *Desulfovibrio desulfuricans* contains numerous HS-producing reactions, whereas *Lactobacillus pasteurii* contains none. Since HS is widespread across microbes, *L. pasteurii* is unusual. This raises the question: how does HS metabolism compare across different *Lactobacillus* species? We can generate a heatmap of reaction presence across reconstruction visualisations using `visualizeNormalizedRxnPresence`. First, download all *Lactobacillus* reconstructions from the Harvard Dataverse. (NB: if the MATLAB Parallel Computing Toolbox is not installed, we recommend downloading only a small subset, e.g. three reconstructions, as sequential computation time can otherwise be excessive.) For simplicity, create a new working directory and place the MicroMap and all pan-*Lactobacillus* reconstructions inside.

```
%First create the reconstructions visualisations
visualizeReconstructionsOnMap('MicroMap.xml', pwd);

%Then create a normalised rxn presence heatmap
visualizeNormalizedRxnPresence(pwd);

%Highlight h2s in the normalised reaction presence heatmap.
[xmlHeatmap, mapHeatmap] = ...
```



```
transformXML2Map('normalisedReactionHeatmap.xml');
mapHeatmap = changeMetColor(mapHeatmap, {'h2s[c]'}, 'LIME');
mapHeatmap = changeNodesArea(mapHeatmap, {'h2s[c]'}, 75, 150);

%Reconstitute the map
transformMap2XML(xmlHeatmap, ...
    mapHeatmap, 'H2SLactobacillusNormRxnHeatmap.xml')
```

From the heatmap, it is evident that HS-involving reactions are not ubiquitous across *Lactobacillus* species; most lack the capacity for HS metabolism (**Figure 4**). This highlights variability in microbial contributions to sulphur cycling in the gut, with *Lactobacillus* playing only a minimal role.




Figure 4: Heatmap of relative reaction abundance for all pan *Lactobacillus* species contained within AGORA2. The thin, dark blue lines indicate a reaction presence of <30%. The HS molecules are highlighted to facilitate inspecting which reactions are attached to HS and therefore engage in HS metabolism.

Visualising flux results from flux balance analysis

The MicroMap enables visualisation of flux results from microbiome metabolic modelling. Flux balance analysis (FBA) can be used to estimate the maximum yields of key metabolites, such as HS in our case study. For a general introduction to FBA, we recommend first consulting the [introductory examples with the E. coli core model](#).

Let's investigate HS secretion in *Desulfovibrio desulfuricans*, a microbe known to play a key role in hydrogen sulphide production.² You can download the *D. desulfuricans* pan-species reconstruction [from the Harvard](#)

[Dataverse](#). Using FBA, we can simulate the maximum secretion of HS. To do this, we set the objective function to the HS exchange reaction and adjust its lower bound to 0, thereby preventing uptake of HS into the cell.

```
%Load the reconstruction of interest
dDesu = readCbModel('panDesulfovibrio_desulfuricans.mat');

%Load and decolorise the MicroMap, if not still loaded from previously
[xmlMicroMap, mapMicroMap] = transformXML2Map('MicroMap.xml');
mapMicroMapUnified = unifyMetabolicMapCD(mapMicroMap);

%Set the objective function to h2s secretion
dDesu = changeObjective(dDesu, 'EX_h2s(e)');
dDesu = changeRxnBounds(dDesu, 'EX_h2s(e)', 0, 'l');
```

In the downloaded reconstruction, exchange reaction bounds are virtually unrestricted. Modelling HS secretion from such an unlimited sulphur supply risks computational artefacts, where results reflect mathematical boundaries rather than biochemical mechanisms. To avoid this, we set uptake of major sulphur-containing compounds to physiologically more realistic levels by constraining the lower bounds of their exchange reactions to around $-20 \text{ mmol}\cdot\text{g}^{-1}\cdot\text{DW}^{-1}$.

```
%Exchange reactions pertaining to most major sulphur-containing compounds
allS = {
    'EX_Lcyst(e)'
    'EX_butso3(e)'
    'EX_cgly(e)'
    'EX_cys_L(e)'
    'EX_ethso3(e)'
    'EX_glycys(e)'
    'EX_glymet(e)'
    'EX_h2s(e)'
    'EX_hexs(e)'
    'EX_isetac(e)'
    'EX_met_D(e)'
    'EX_met_L(e)'
    'EX_metala(e)'
    'EX_metsox_R_L(e)'
    'EX_metsox_S_L(e)'
    'EX_mops(e)'
    'EX_mso3(e)'
    'EX_so4(e)'
    'EX_sulfac(e)'
    'EX_taur(e)'
    'EX_tchola(e)'
    'EX_tdchola(e)'
    'EX_tdechola(e)'
    'EX_tsul(e)'
};

%Set the uptake rate for all major sulfur-containin compounds
```



```
dDesu = changeRxnBounds(dDesu, allS, -20, 'l');
```

We can then run FBA to estimate the maximum HS secretion under these conditions. The associated flux distribution can be visualised with the *addFluxFBA* function.

```
%Solve FBA and visualise the results
allSFBA = optimizeCbModel(dDesu);

mapAllSFBA = addFluxFBA(mapMicroMapUnified, dDesu, allSFBA, 'DARKVIOLET');
transformMap2XML(xmlMicroMap, mapAllSFBA, 'H2SMicroMapAllSFBA.xml');
```

We can compare results by running FBA under different conditions. From our earlier MicroMap inspection, HS production is linked to both inorganic sulphur compounds and sulphur-containing amino acids such as cysteine and methionine. Accordingly, we shall run FBA with uptake permitted primarily from inorganic compounds, and then from cysteine and methionine.

```
%Define list of inorganic sulphur-containing compounds
inorgs = {
    'EX_tsul(e)'
    'EX_so4(e)'
}

%Define list of cysteine and methionine compounds
cysMets = {
    'EX_cgly(e)'
    'EX_cys_L(e)'
    'EX_glycys(e)'
    'EX_glymet(e)'
    'EX_met_D(e)'
    'EX_met_L(e)'
    'EX_metala(e)'
}

%Set all uptake bounds for the major sulphur-containing compounds to 0
dDesu = changeRxnBounds(dDesu, allS, 0, 'l');

%Set the uptake rate for inorganic sulphur compounds
dDesu = changeRxnBounds(dDesu, inorgs, -20, 'l');

%Solve FBA and visualise the results
inorgsFBA = optimizeCbModel(dDesu);

mapInorgsFBA = addFluxFBA(mapMicroMapUnified, dDesu, inorgsFBA,
'DARKVIOLET');
transformMap2XML(xmlMicroMap, mapInorgsFBA, 'H2SMicroMapInorgsFBA.xml');

%Set the uptake rate for cysteine and methionine sulphur compounds
dDesu = changeRxnBounds(dDesu, inorgs, 0, 'l');
```

```
dDesu = changeRxnBounds(dDesu, cysMets, -20, '1');

%Solve FBA and visualise the results
cysMetsFBA = optimizeCbModel(dDesu);

mapCysMetsFBA = addFluxFBA(mapMicroMapUnified, dDesu, cysMetsFBA,
'DARKVIOLET');
transformMap2XML(xmlMicroMap, mapCysMetsFBA, 'H2SMicroMapCysMetsFBA.xml');
```

Inspection of the flux vectors shows that HS secretion is higher on a cysteine–methionine medium (100 mmol·g⁻¹·DW⁻¹) than on an inorganic sulphur donor medium (60 mmol·g⁻¹·DW⁻¹). Other sulphur compounds also contribute, as indicated by the full donor set (200 mmol·g⁻¹·DW⁻¹) (**Figure 5**). Remember that, while the FBA-derived flux result is unique, the underlying flux distribution may not be, since the system of linear equations is underdetermined. Nonetheless, these distributions remain valuable for interpreting the systems context of the solution.






Figure 5: Comparison of FBA results obtained under different conditions. The allSFBA result allowed uptake from all major sulphur sources, cysMetsFBA allowed from cysteine and methionine, and inorgsFBA allowed from inorganic sulphur sources.

Note that all the solutions involve many different reactions across the MicroMap, well beyond sulphur metabolism. Comparing the results from the inorganics and cysMets FBA, we see that the cysMets solution does not carry flux in the region of the 'Sulfur Metabolism' biochemical subsystem, which highlights several reactions deriving from inorganic sulphur compounds (**Figure 6**).



Cysteine & methionine



Inorganics

Figure 6: Comparison of flux distributions from FBA for maximising H₂S secretion derived from using cysteine and methionine versus inorganic sulphur sources. No flux is observed in the ‘Sulfur Metabolism’ subsystem from a system of cysteine and methionine sulphur sources, whereas flux is observed from a system of inorganic sulphur sources.

Another interesting area is around the methionine subsystem, where the cysMets solution carries flux in more reactions. Upon closer inspection, one sees that the inorganics solution carries flux in reactions synthesising cysteine (**Figure 7**). The cysMets solution, does not carry flux in that same reaction to synthesise cysteine, but conversely, carries flux in reactions synthesising thiosulphate (**Figure 7**). This flux distribution highlights the interconnectedness and redundance for some metabolite syntheses.





Figure 7: Comparison of flux distributions from FBA results for maximising H₂S secretion derived from using cysteine and methionine versus inorganic sulphur sources. Cysteine can be produced from a system of inorganic sulphur sources, whereas thiosulphate can be produced from a system of cysteine and methionine sulphur sources.

In this case study, we have used the MicroMap to follow hydrogen sulphide metabolism step by step: first through manual inspection in CellDesigner, and then programmatically via the COBRA Toolbox. We learned that some microbes, such as *Desulfovibrio desulfuricans*, possess extensive capacity for HS production, while others, such as *Lactobacillus* species, contribute little. The MicroMap can be flexibly integrated into a wide range of workflows, and we encourage readers to adapt and extend these approaches to their own research questions.

References

1. Heirendt L, Arreckx S, Pfau T, Mendoza SN, Richelle A, Heinken A, et al. Creation and analysis of biochemical constraint-based models using the COBRA Toolbox v.3.0. *Nat Protoc.* 2019;14(3):639-702.
2. Mutuyemungu E, Singh M, Liu S, Rose DJ. Intestinal gas production by the gut microbiota: A review. *Journal of Functional Foods.* 2023;100:105367.
3. Suarez FL, Springfield J, Levitt MD. Identification of gases responsible for the odour of human flatus and evaluation of a device purported to reduce this odour. *Gut.* 1998;43(1):100-4.
4. Kushkevych I, Dordevi D, Kollar P, Vít zová M, Drago L. Hydrogen Sulfide as a Toxic Product in the Small–Large Intestine Axis and its Role in IBD Development. *Journal of Clinical Medicine.* 2019;8(7):1054.